

TEM Tissue Processing

Prior to starting an Electron Microscopy Core project, investigator must meet with Senior Scientist Anza Darehshouri to discuss project goals, appropriate processing protocols, microscopes and request fixation solutions for preserving tissues.

The facility is unable to accept tissue which has not been fixed. Please mince tissue to small size (approx. 1 mm square or smaller). **Ensure that tissue does not dry out during this time period. Artifacts will be evident in the TEM if tissue is allowed to dry.**

Routine TEM tissue processing usually takes three days from process start date. Additionally, ultramicrotomy and grid staining must also be done before the specimen is ready for the scope. **Thus, the minimum turn around time for standard processing is no less than one week from date of submission.**

Day 1

1. Primary Fixation 2.5% Glutaraldehyde in 0.1M Cacodylate buffer, pH 7.4 @RT for 1.5 hours, or overnight 4°C
2. Rinse 5X 0.1M Cacodylate buffer pH7.4 @ RT 5-8 minutes each
3. Secondary Fixation 1-2X 1% Osmium Tetroxide with 0.8% K₃Fe (CN)₆ in cacodylate buffer for 1.5 hours RT in dark
4. Rinse 5X Double distilled water @RT 5-8 minutes each (last wash overnight @4°C)

Day2

1. Pre-Stain 4% Uranyl Acetate in 50% ETOH @ RT 2 hours in dark
2. Dehydration 50% Ethanol in DH₂O 1X @RT for 8 minutes
70% Ethanol in DH₂O 2X @RT for 8 minutes
85% Ethanol in DH₂O 1X @RT for 8 minutes
95% Ethanol in DH₂O 2X @RT for 8 minutes
100% in DH₂O 4X @RT for 8 minutes (Open new 100% ETOH bottle if more than 1 week old)
3. Transition Propylene Oxide (PO) 1X @RT for 10 minutes
4. Infiltration PO:EPON (2:1) 1X @RT for 2+ hours on rotator
PO:EPON (1:2) 1X @RT for 2+ hours on rotator then overnight

Day3

1. Infiltration 100% EPON 2X @RT for 2+ hours on rotator
2. Embed 100% EPON Add paper label with sample ID to embedding mold
3. Polymerization Place embedding molds in 70°C oven overnight

Date:
Researcher:
P.I.:
Specimen: