

IMPLANT RUPTURE: Silicone

Anatomy: Implants
 Sub-Anatomy: Breast

- Exams
 - Routine

Coil: 7 Channel Breast Coil

SEQUENCE - BASICS				VARIABLE PARAMETERS													
PLANE	SEQ	Slice / Gap (mm)	Misc. / Comment	MTX	FOV (cm)	% RFO V	Scan %	Pixel size (mm)	TR	TE	FA/TI	Nex Avg Acq	NS	ETL Turbo Factor	Phase Encode	Scan TIME (min)	Pixel Shift BW-kHz
ROUTINE																	
AX	T1 VISTA 3D	200	FOV PATIENT DEPEND						2000	204							
AX	T1W 3D	200	FOV PATIENT DEPEND						400	15							
AX	Silicone Bright	200	FOV PATIENT DEPEND						13519/ 2700	70							
AX	Silicone Suppression	200	FOV PATIENT DEPEND						7.3	3.6							
SAG	T2FS LEFT	30	FOV PATIENT DEPEND						3708	70							
SAG	T2FS RIGHT	30	FOV PATIENT DEPEND						3708	70							
COR	T2 W/TSE	50	FOV PATIENT DEPEND						2867	140							

TIPS: