

UT Southwestern Department of Radiology

Protocol Name: Pancreas 2 Phase plus Pelvis

Orderable Name: CT ABDOMEN/PELVIS W IV CONTRAST

Adult only

Protocol Epic Button: Panc 2 + Pelvis

CTDIvol < 60 mGy
DLP < 10000 mGycm

Indications: Pancreatic mass (known or suspected, initial evaluation, resectability, or follow-up). Use for all pancreatic masses (adenocarcinoma or neuroendocrine) and pancreatic cyst evaluation. For evaluation of pancreatitis, use abdomen 1 or AP1

Active Protocol

Acquisitions: 2

<p>Oral Contrast: Water</p> <p>H2O Variant:</p> <p>Oral Effervescent Granules on Table? <input type="checkbox"/></p> <p>Regime: Notes:</p> <p>12 oz given 30 min prior to exam 12 oz given on the table</p>	<p>IV Contrast: Link to Contrast Information</p> <p>Rate (mL/sec): P3T Volume (mL): P3T</p> <p>IV Access: Power injection: 20g or larger in large vein (prefer AC fossa or forearm)</p> <p>Notes:</p>	<p>Other Contrast: None</p> <p>Airway: Full inspiration</p> <hr/> <p>Other Notes: Consult body habitus kVp selection chart.</p>
--	--	--

Last Change: 1/19/2022

Last Review: 1/19/2022

Links: [kVp Body Chart](#)

[General Statements](#)

Special Instructions		
Series Name	1 Late Arterial	2 Portal Venous
Phase Timing	Bolus tracking	Bolus tracking
Acquisition Protocol	Abdominal Contrast 3mm	Abdominal Contrast 3mm
Coverage	Diaphragm to iliac crests	Diaphragm thru pubis symph
FOV	Whole body	Whole body
Algorithm	Soft Tissue	Soft Tissue
Primary Axial Recon	3 mm	3 mm
Other Recons	2 mm coronal and sagittal	2 mm coronal and sagittal
MIP Recons		
*Dual Energy Philips	mono E 50, SBI	mono E 50, SBI
*Dual Energy Siemens	VNC, Iodine map, mono E 50	VNC, Iodine map, mono E 50

* When dual energy or spectral CT is used

