

Anatomy: **ABDOMEN**
 Sub-Anatomy: **Non-contrast Abdomen non BH**

- Exams
 - Routine Coil: **SENSE XL Torso**

SEQUENCE - BASICS				VARIABLE PARAMETERS													
PLANE	SEQ	Slice / Gap (mm)	Misc / Comment	MTX	FOV (cm)	% RFOV	Scan %	Pixel size (mm)	TR	TE	FA/TI	Nex Avg Acq	NS	ETL Turbo Factor	Phase Encode	Scan TIME (min)	Pixel Shift BW-kHz
ROUTINE																	
COR	2D T2 SSFSE	5/1	Resp Trig														
AX	2D T2 SSFSE	5/1	Resp Trig														
AX	2D T2 FS SSFSE	5/1	Resp Trig														
AX	2D T1 TFE IP	5/1	Resp Trig														
AX	2D T1 TFE OP	5/1	Resp Trig														
AX	2D DWI FS	7/1	b=0,50,400,800														
↓ OPTIONAL ↓																	

ORDERABLE: MR ABDOMEN WO CONTRAST.

TIPS:

